

Antibodies Product Datasheet

Product Name:	EphA1 Monoclonal Antibody
Catalog Number:	N0219
Reactivity:	Human
Applications:	WB,IHC-p,IF(paraffin section),ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of EphA1 expressed in E. Coli.
Concentration:	1 mg/ml
Gene Name:	EPHA1
Protein Name:	Ephrin type-A receptor 1
Other Name:	EPHA1; EPH; EPHT; EPHT1; Ephrin type-A receptor 1; hEpha1; EPH tyrosine kinase; EPH tyrosine kinase 1; Erythropoietin-producing hepatoma receptor; Tyrosine-protein kinase receptor EPH
Human Swiss Prot No.:	P21709
Formulation:	Purified antibody in PBS containing 0.03% sodium azide.
Specificity:	EphA1 Monoclonal Antibody detects endogenous levels of EphA1 protein.
Storage and Stability:	-20°C/1 year