

Antibodies Product Datasheet

Product Name:	EphA2 Monoclonal Antibody
Catalog Number:	N0221
Reactivity:	Human
Applications:	WB,IHC-p,IF(paraffin section),ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of EphA2 expressed in E. Coli.
Concentration:	
Gene Name:	EPHA2
Protein Name:	Ephrin type-A receptor 2
Other Name:	EPHA2; ECK; Ephrin type-A receptor 2; Epithelial cell kinase; Tyrosine-protein kinase receptor ECK
Human Swiss Prot No.:	P29317
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	EphA2 Monoclonal Antibody detects endogenous levels of EphA2 protein.
Storage and Stability:	-20°C/1 year