

Antibodies Product Datasheet

Product Name:	EphA4 Monoclonal Antibody
Catalog Number:	N0223
Reactivity:	Human
Applications:	WB,ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of EphA4 (aa777-986) expressed in E. Coli.
Concentration:	
Gene Name:	EPHA4
Protein Name:	Ephrin type-A receptor 4
Other Name:	EPHA4; HEK8; SEK; TYRO1; Ephrin type-A receptor 4; EPH-like kinase 8; EK8; hEK8; Tyrosine-protein kinase TYRO1; Tyrosine-protein kinase receptor SEK
Human Swiss Prot No.:	P54764
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	EphA4 Monoclonal Antibody detects endogenous levels of EphA4 protein.
Storage and Stability:	-20°C/1 year