

Antibodies Product Datasheet

Product Name:	EphA6 Monoclonal Antibody
Catalog Number:	N0225
Reactivity:	Human
Applications:	WB,ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of EphA6 (aa695-795) expressed in E. Coli.
Concentration:	
Gene Name:	EPHA6
Protein Name:	Ephrin type-A receptor 6
Other Name:	EPHA6; EHK2; HEK12; Ephrin type-A receptor 6; EPH homology kinase 2; EHK-2; EPH-like kinase 12; EK12
Human Swiss Prot No.:	Q9UF33
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	EphA6 Monoclonal Antibody detects endogenous levels of EphA6 protein.
Storage and Stability:	-20°C/1 year