

Antibodies Product Datasheet

Product Name:	EphA7 Monoclonal Antibody
Catalog Number:	N0226
Reactivity:	Human
Applications:	WB,IHC-p,IF(paraffin section),ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of EphA7 (aa27-210) expressed in E. Coli.
Concentration:	
Gene Name:	EPHA7
Protein Name:	Ephrin type-A receptor 7
Other Name:	EPHA7; EHK3; HEK11; Ephrin type-A receptor 7; EPH homology kinase 3; EHK-3; EPH-like kinase 11; EK11; hEK11
Human Swiss Prot No.:	Q15375
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	EphA7 Monoclonal Antibody detects endogenous levels of EphA7 protein.
Storage and Stability:	-20°C/1 year