

Antibodies Product Datasheet

Product Name:	EphA8 Monoclonal Antibody
Catalog Number:	N0227
Reactivity:	Human
Applications:	WB,ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of EphA8 (aa70-150) expressed in E. Coli.
Concentration:	
Gene Name:	EPHA8
Protein Name:	Ephrin type-A receptor 8
Other Name:	EPHA8; EEK; HEK3; KIAA1459; Ephrin type-A receptor 8; EPH- and ELK-related kinase; EPH-like kinase 3; EK3; hEK3; Tyrosine-protein kinase receptor EEK
Human Swiss Prot No.:	P29322
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	EphA8 Monoclonal Antibody detects endogenous levels of EphA8 protein.
Storage and Stability:	-20°C/1 year