

Antibodies Product Datasheet

Product Name:	EphB2 Monoclonal Antibody
Catalog Number:	N0229
Reactivity:	Human
Applications:	WB,IF/ICC,ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of EphB2 (aa17-200) expressed in E. Coli.
Concentration:	
Gene Name:	EPHB2
Protein Name:	Ephrin type-B receptor 2
Other Name:	EPHB2; DRT; EPHT3; EPTH3; ERK; HEK5; TYRO5; Ephrin type-B receptor 2; Developmentally-regulated Eph-related tyrosine kinase; ELK-related tyrosine kinase; EPH tyrosine kinase 3; EPH-like kinase 5; EK5; hEK5; Renal carcinoma antigen NY-REN-47
Human Swiss Prot No.:	P29323
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	EphB2 Monoclonal Antibody detects endogenous levels of EphB2 protein.
Storage and Stability:	-20°C/1 year