

Antibodies Product Datasheet

Product Name:	EphB3 Monoclonal Antibody
Catalog Number:	N0230
Reactivity:	Human
Applications:	WB,IHC-p,IF(paraffin section),ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of EphB3 (aa39-212) expressed in E. Coli.
Concentration:	
Gene Name:	EPHB3
Protein Name:	Ephrin type-B receptor 3
Other Name:	EPHB3; ETK2; HEK2; TYRO6; Ephrin type-B receptor 3; EPH-like tyrosine kinase 2; EPH-like kinase 2; Embryonic kinase 2; EK2; hEK2; Tyrosine-protein kinase TYRO6
Human Swiss Prot No.:	P54753
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	EphB3 Monoclonal Antibody detects endogenous levels of EphB3 protein.
Storage and Stability:	-20°C/1 year