

Antibodies Product Datasheet

Product Name:	ERK 3 Monoclonal Antibody
Catalog Number:	N0245
Reactivity:	Human
Applications:	WB,FCM,ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Flow cytometry: 1/200 - 1/400. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human ERK 3 expressed in E. Coli.
Concentration:	
Gene Name:	MAPK6
Protein Name:	Mitogen-activated protein kinase 6
Other Name:	MAPK6; ERK3; PRKM6; Mitogen-activated protein kinase 6; MAP kinase 6; MAPK 6; Extracellular signal-regulated kinase 3; ERK-3; MAP kinase isoform p97; p97-MAPK
Human Swiss Prot No.:	Q16659
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	ERK 3 Monoclonal Antibody detects endogenous levels of ERK 3 protein.
Storage and Stability:	-20°C/1 year