

Antibodies Product Datasheet

Product Name:	ER α Monoclonal Antibody
Catalog Number:	N0251
Reactivity:	Human
Applications:	WB,IHC-p,IF/ICCELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/200 - 1/1000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human ER α expressed in E. Coli.
Concentration:	
Gene Name:	ESR1
Protein Name:	Estrogen receptor
Other Name:	ESR1; ESR; NR3A1; Estrogen receptor; ER; ER-alpha; Estradiol receptor; Nuclear receptor subfamily 3 group A member 1
Human Swiss Prot No.:	P03372
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	ER α Monoclonal Antibody detects endogenous levels of ER α protein.
Storage and Stability:	-20°C/1 year