

Antibodies Product Datasheet

Product Name:	ESET Monoclonal Antibody
Catalog Number:	N0253
Reactivity:	Human, Mouse, Monkey
Applications:	WB, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human ESET expressed in E. Coli.
Concentration:	
Gene Name:	SETDB1
Protein Name:	Histone-lysine N-methyltransferase SETDB1
Other Name:	SETDB1; KIAA0067; KMT1E; Histone-lysine N-methyltransferase SETDB1; ERG-associated protein with SET domain; ESET; Histone H3- K9 methyltransferase 4; H3-K9-HMTase 4; Lysine N-methyltransferase 1E; SET domain bifurcated 1
Human Swiss Prot No.:	Q15047
Formulation:	Ascitic fluid containing 0.03% sodium azide, 0.5% BSA, 50% glycerol.
Specificity:	ESET Monoclonal Antibody detects endogenous levels of ESET protein.
Storage and Stability:	-20°C/1 year