

Antibodies Product Datasheet

Product Name:	Factor VIII Monoclonal Antibody
Catalog Number:	N0256
Reactivity:	Human
Applications:	WB,ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of Factor VIII expressed in E. Coli.
Concentration:	
Gene Name:	F8
Protein Name:	Coagulation factor VIII
Other Name:	F8; F8C; Coagulation factor VIII; Antihemophilic factor; AHF; Procoagulant component
Human Swiss Prot No.:	P00451
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	Factor VIII Monoclonal Antibody detects endogenous levels of Factor VIII protein.
Storage and Stability:	-20°C/1 year