

Antibodies Product Datasheet

Product Name:	FAK Monoclonal Antibody
Catalog Number:	N0258
Reactivity:	Human, Mouse
Applications:	WB, IHC-p, IF/ICELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/200 - 1/1000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human FAK expressed in E. Coli.
Concentration:	
Gene Name:	PTK2
Protein Name:	Focal adhesion kinase 1
Other Name:	PTK2; FAK; FAK1; Focal adhesion kinase 1; FADK 1; Focal adhesion kinase-related nonkinase; FRNK; Protein phosphatase 1 regulatory subunit 71; PPP1R71; Protein-tyrosine kinase 2; p125FAK; pp125FAK
Human Swiss Prot No.:	Q05397
Formulation:	Ascitic fluid containing 0.03% sodium azide, 0.5% BSA, 50% glycerol.
Specificity:	FAK Monoclonal Antibody detects endogenous levels of FAK protein.
Storage and Stability:	-20°C/1 year