

Antibodies Product Datasheet

Product Name:	Fer Monoclonal Antibody
Catalog Number:	N0261
Reactivity:	Human, Mouse
Applications:	WB, IHC-p, IF/ICELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/200 - 1/1000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human Fer expressed in E. Coli.
Concentration:	
Gene Name:	FER
Protein Name:	Tyrosine-protein kinase Fer
Other Name:	FER; TYK3; Tyrosine-protein kinase Fer; Feline encephalitis virus-related kinase FER; Fujinami poultry sarcoma/Feline sarcoma-related protein Fer; Proto-oncogene c-Fer; Tyrosine kinase 3; p94-Fer
Human Swiss Prot No.:	P16591
Formulation:	Ascitic fluid containing 0.03% sodium azide, 0.5% BSA, 50% glycerol.
Specificity:	Fer Monoclonal Antibody detects endogenous levels of Fer protein.
Storage and Stability:	-20°C/1 year