

Antibodies Product Datasheet

Product Name:	Fes Monoclonal Antibody
Catalog Number:	N0262
Reactivity:	Human
Applications:	WB,IHC-p,IF(paraffin section),ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of Fer expressed in E. Coli.
Concentration:	
Gene Name:	FES
Protein Name:	Tyrosine-protein kinase Fes/Fps
Other Name:	FES; FPS; Tyrosine-protein kinase Fes/Fps; Feline sarcoma/Fujinami avian sarcoma oncogene homolog; Proto-oncogene c-Fes; Proto- oncogene c-Fps; p93c-fes
Human Swiss Prot No.:	P07332
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	Fes Monoclonal Antibody detects endogenous levels of Fes protein.
Storage and Stability:	-20°C/1 year