

Antibodies Product Datasheet

Product Name:	Fibulin-5 Monoclonal Antibody
Catalog Number:	N0271
Reactivity:	Human
Applications:	WB,IHC-p,IF(paraffin section),ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of Fibulin-5 expressed in E. Coli.
Concentration:	
Gene Name:	FBLN5
Protein Name:	Fibulin-5
Other Name:	FBLN5; DANCE; Fibulin-5; FIBL-5; Developmental arteries and neural crest EGF-like protein; Dance; Urine p50 protein; UP50
Human Swiss Prot No.:	Q9UBX5
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	Fibulin-5 Monoclonal Antibody detects endogenous levels of Fibulin-5 protein.
Storage and Stability:	-20°C/1 year