

Antibodies Product Datasheet

Product Name:	Flt-1 Monoclonal Antibody
Catalog Number:	N0277
Reactivity:	Human
Applications:	WB,ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant extracellular fragment of human Flt-1 fused with hlgGfc tag expressed in HEK293 cells.
Concentration:	
Gene Name:	FLT1
Protein Name:	Vascular endothelial growth factor receptor 1
Other Name:	FLT1; FLT; FRT; VEGFR1; Vascular endothelial growth factor receptor 1; VEGFR-1; Fms-like tyrosine kinase 1; FLT-1; Tyrosine-protein kinase FRT; Tyrosine-protein kinase receptor FLT; FLT; Vascular permeability factor receptor
Human Swiss Prot No.:	P17948
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	Flt-1 Monoclonal Antibody detects endogenous levels of Flt-1 protein.
Storage and Stability:	-20°C/1 year