

Antibodies Product Datasheet

Product Name:	FoxP2 Monoclonal Antibody
Catalog Number:	N0283
Reactivity:	Human
Applications:	WB,ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human FoxP2 expressed in E. Coli.
Concentration:	
Gene Name:	FOXP2
Protein Name:	Forkhead box protein P2
Other Name:	FOXP2; CAGH44; TNRC10; Forkhead box protein P2; CAG repeat protein 44; Trinucleotide repeat-containing gene 10 protein
Human Swiss Prot No.:	O15409
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	FoxP2 Monoclonal Antibody detects endogenous levels of FoxP2 protein.
Storage and Stability:	-20°C/1 year