

Antibodies Product Datasheet

Product Name:	FoxP3 Monoclonal Antibody
Catalog Number:	N0285
Reactivity:	Human, Mouse
Applications:	WB, IHC-p, IF/ICELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/200 - 1/1000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human FoxP3 expressed in E. Coli.
Concentration:	1 mg/ml
Gene Name:	FOXP3
Protein Name:	Forkhead box protein P3
Other Name:	FOXP3; IPEX; JM2; Forkhead box protein P3; Scurfin
Human Swiss Prot No.:	Q9BZS1
Formulation:	Purified antibody in PBS containing 0.03% sodium azide.
Specificity:	FoxP3 Monoclonal Antibody detects endogenous levels of FoxP3 protein.
Storage and Stability:	-20°C/1 year