

Antibodies Product Datasheet

Product Name:	Fucokinase Monoclonal Antibody
Catalog Number:	N0287
Reactivity:	Human, Mouse, Rat, Monkey
Applications:	WB, FCM, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Flow cytometry: 1/200 - 1/400. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human Fucokinase expressed in E. Coli.
Concentration:	
Gene Name:	FUK
Protein Name:	L-fucose kinase
Other Name:	FUK; L-fucose kinase; Fucokinase
Human Swiss Prot No.:	Q8N0W3
Formulation:	Ascitic fluid containing 0.03% sodium azide, 0.5% BSA, 50% glycerol.
Specificity:	Fucokinase Monoclonal Antibody detects endogenous levels of Fucokinase protein.
Storage and Stability:	-20°C/1 year