

Antibodies Product Datasheet

Product Name:	Fyn Monoclonal Antibody
Catalog Number:	N0289
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human Fyn expressed in E. Coli.
Concentration:	
Gene Name:	FYN
Protein Name:	Tyrosine-protein kinase Fyn
Other Name:	FYN; Tyrosine-protein kinase Fyn; Proto-oncogene Syn; Proto-oncogene c-Fyn; Src-like kinase; SLK; p59-Fyn
Human Swiss Prot No.:	P06241
Formulation:	Ascitic fluid containing 0.03% sodium azide, 0.5% BSA, 50% glycerol.
Specificity:	Fyn Monoclonal Antibody detects endogenous levels of Fyn protein.
Storage and Stability:	-20°C/1 year