

Antibodies Product Datasheet

Product Name:	GABP- α Monoclonal Antibody
Catalog Number:	N0292
Reactivity:	Human, Mouse
Applications:	WB, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human GABP- α (aa120-190) expressed in E. Coli.
Concentration:	
Gene Name:	GABPA
Protein Name:	GA-binding protein alpha chain
Other Name:	GABPA; E4TF1A; GA-binding protein alpha chain; GABP subunit alpha; Nuclear respiratory factor 2 subunit alpha; Transcription factor E4TF1-60
Human Swiss Prot No.:	Q06546
Formulation:	Ascitic fluid containing 0.03% sodium azide, 0.5% BSA, 50% glycerol.
Specificity:	GABP- α Monoclonal Antibody detects endogenous levels of GABP- α protein.
Storage and Stability:	-20°C/1 year