

Antibodies Product Datasheet

Product Name:	GATA-1 Monoclonal Antibody
Catalog Number:	N0295
Reactivity:	Human
Applications:	WB,IHC-p,IF/ICCELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/200 - 1/1000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human GATA-1 expressed in E. Coli.
Concentration:	
Gene Name:	GATA1
Protein Name:	Erythroid transcription factor
Other Name:	GATA1; ERYF1; GF1; Erythroid transcription factor; Eryf1; GATA-binding factor 1; GATA-1; GF-1; NF-E1 DNA-binding protein
Human Swiss Prot No.:	P15976
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	GATA-1 Monoclonal Antibody detects endogenous levels of GATA-1 protein.
Storage and Stability:	-20°C/1 year