

Antibodies Product Datasheet

Product Name:	GATA-4 Monoclonal Antibody
Catalog Number:	N0298
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human GATA-4 expressed in E. Coli.
Concentration:	
Gene Name:	GATA4
Protein Name:	Transcription factor GATA-4
Other Name:	GATA4; Transcription factor GATA-4; GATA-binding factor 4
Human Swiss Prot No.:	P43694
Formulation:	Ascitic fluid containing 0.03% sodium azide, 0.5% BSA, 50% glycerol.
Specificity:	GATA-4 Monoclonal Antibody detects endogenous levels of GATA-4 protein.
Storage and Stability:	-20°C/1 year