

Antibodies Product Datasheet

Product Name:	GATA-5 Monoclonal Antibody
Catalog Number:	N0299
Reactivity:	Human
Applications:	WB,ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human GATA-5 expressed in E. Coli.
Concentration:	
Gene Name:	GATA5
Protein Name:	Transcription factor GATA-5
Other Name:	GATA5; Transcription factor GATA-5; GATA-binding factor 5
Human Swiss Prot No.:	Q9BWX5
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	GATA-5 Monoclonal Antibody detects endogenous levels of GATA-5 protein.
Storage and Stability:	-20°C/1 year