

Antibodies Product Datasheet

Product Name:	Gfi-1 Monoclonal Antibody
Catalog Number:	N0304
Reactivity:	Human
Applications:	WB,ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human Gfi-1 expressed in E. Coli.
Concentration:	
Gene Name:	GFI1
Protein Name:	Zinc finger protein Gfi-1
Other Name:	GFI1; ZNF163; Zinc finger protein Gfi-1; Growth factor independent protein 1; Zinc finger protein 163
Human Swiss Prot No.:	Q99684
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	Gfi-1 Monoclonal Antibody detects endogenous levels of Gfi-1 protein.
Storage and Stability:	-20°C/1 year