

Antibodies Product Datasheet

Product Name:	GKLF Monoclonal Antibody
Catalog Number:	N0306
Reactivity:	Human
Applications:	WB,IHC-p,IF(paraffin section),ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human GKLF expressed in E. Coli.
Concentration:	
Gene Name:	MS4A1
Protein Name:	B-lymphocyte antigen CD20
Other Name:	MS4A1; CD20; B-lymphocyte antigen CD20; B-lymphocyte surface antigen B1; Bp35; Leukocyte surface antigen Leu-16; Membrane-spanning 4-domains subfamily A member 1; CD antigen CD20
Human Swiss Prot No.:	P11836
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	GKLF Monoclonal Antibody detects endogenous levels of GKLF protein.
Storage and Stability:	-20°C/1 year