

Antibodies Product Datasheet

Product Name:	GluR-3 Monoclonal Antibody
Catalog Number:	N0308
Reactivity:	Human
Applications:	WB,IHC-p,IF(paraffin section),ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human GluR-3 expressed in E. Coli.
Concentration:	
Gene Name:	GRIA3
Protein Name:	Glutamate receptor 3
Other Name:	GRIA3; GLUR3; GLURC; Glutamate receptor 3; GluR-3; AMPA-selective glutamate receptor 3; GluR-C; GluR-K3; Glutamate receptor ionotropic; AMPA 3; GluA3
Human Swiss Prot No.:	P42263
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	GluR-3 Monoclonal Antibody detects endogenous levels of GluR-3 protein.
Storage and Stability:	-20°C/1 year