

Antibodies Product Datasheet

Product Name:	GPI Monoclonal Antibody
Catalog Number:	N0311
Reactivity:	Human
Applications:	WB,IHC-p,IF/ICCELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/200 - 1/1000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human GPI expressed in E. Coli.
Concentration:	
Gene Name:	GPI
Protein Name:	Glucose-6-phosphate isomerase
Other Name:	GPI; Glucose-6-phosphate isomerase; GPI; Autocrine motility factor; AMF; Neuroleukin; NLK; Phosphoglucose isomerase; PGI; Phosphohexose isomerase; PHI; Sperm antigen 36; SA-36
Human Swiss Prot No.:	P06744
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	GPI Monoclonal Antibody detects endogenous levels of GPI protein.
Storage and Stability:	-20°C/1 year