

Antibodies Product Datasheet

Product Name:	GRK 2 Monoclonal Antibody
Catalog Number:	N0313
Reactivity:	Human, Mouse, Rat, Monkey
Applications:	WB, IHC-p, IF/ICELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/200 - 1/1000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human GRK 2 expressed in E. Coli.
Concentration:	
Gene Name:	ADRBK1
Protein Name:	Beta-adrenergic receptor kinase 1
Other Name:	ADRBK1; BARK; BARK1; GRK2; Beta-adrenergic receptor kinase 1; Beta-ARK-1; G-protein coupled receptor kinase 2
Human Swiss Prot No.:	P25098
Formulation:	Ascitic fluid containing 0.03% sodium azide, 0.5% BSA, 50% glycerol.
Specificity:	GRK 2 Monoclonal Antibody detects endogenous levels of GRK 2 protein.
Storage and Stability:	-20°C/1 year