

Antibodies Product Datasheet

Product Name:	GSK3 α Monoclonal Antibody
Catalog Number:	N0314
Reactivity:	Human
Applications:	WB,IHC-p,IF(paraffin section),ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of GSK3 α expressed in E. Coli.
Concentration:	
Gene Name:	GSK3A
Protein Name:	Glycogen synthase kinase-3 alpha
Other Name:	GSK3A; Glycogen synthase kinase-3 alpha; GSK-3 alpha; Serine/threonine-protein kinase GSK3A
Human Swiss Prot No.:	P49840
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	GSK3 α Monoclonal Antibody detects endogenous levels of GSK3 α protein.
Storage and Stability:	-20°C/1 year