

Antibodies Product Datasheet

Product Name:	GSTP1 Monoclonal Antibody
Catalog Number:	N0319
Reactivity:	Human
Applications:	WB,IHC-p,IF/ICCFCM,ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/200 - 1/1000. Immunofluorescence: 1/200 - 1/1000. Flow cytometry: 1/200 - 1/400. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human GSTP1 expressed in E. Coli.
Concentration:	
Gene Name:	GSTP1
Protein Name:	Glutathione S-transferase P
Other Name:	GSTP1; FAES3; GST3; Glutathione S-transferase P; GST class-pi; GSTP1-1
Human Swiss Prot No.:	P09211
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	GSTP1 Monoclonal Antibody detects endogenous levels of GSTP1 protein.
Storage and Stability:	-20°C/1 year