

Antibodies Product Datasheet

Product Name:	GTBP Monoclonal Antibody
Catalog Number:	N0320
Reactivity:	Human
Applications:	WB,ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of GTBP expressed in E. Coli.
Concentration:	
Gene Name:	MSH6
Protein Name:	DNA mismatch repair protein Msh6
Other Name:	MSH6; GTBP; DNA mismatch repair protein Msh6; hMSH6; G/T mismatch-binding protein; GTBP; GTMBP; MutS-alpha 160 kDa subunit; p160
Human Swiss Prot No.:	P52701
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	GTBP Monoclonal Antibody detects endogenous levels of GTBP protein.
Storage and Stability:	-20°C/1 year