

Antibodies Product Datasheet

Product Name:	HAND1 Monoclonal Antibody
Catalog Number:	N0322
Reactivity:	Human
Applications:	WB,ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of HAND1 (aa90-190) expressed in E. Coli.
Concentration:	
Gene Name:	HAND1
Protein Name:	Heart- and neural crest derivatives-expressed protein 1
Other Name:	HAND1; BHLHA27; EHAND; Heart- and neural crest derivatives-expressed protein 1; Class A basic helix-loop-helix protein 27; bHLHa27; Extraembryonic tissues; heart, autonomic nervous system and neural crest derivatives-expressed protein 1; eH
Human Swiss Prot No.:	O96004
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	HAND1 Monoclonal Antibody detects endogenous levels of HAND1 protein.
Storage and Stability:	-20°C/1 year