

Antibodies Product Datasheet

Product Name:	Hck Monoclonal Antibody
Catalog Number:	N0325
Reactivity:	Human
Applications:	WB,IHC-p,IF(paraffin section),ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of Hck expressed in E. Coli.
Concentration:	
Gene Name:	HCK
Protein Name:	Tyrosine-protein kinase HCK
Other Name:	HCK; Tyrosine-protein kinase HCK; Hematopoietic cell kinase; Hemopoietic cell kinase; p59-HCK/p60-HCK; p59Hck; p61Hck
Human Swiss Prot No.:	P08631
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	Hck Monoclonal Antibody detects endogenous levels of Hck protein.
Storage and Stability:	-20°C/1 year