

Antibodies Product Datasheet

Product Name:	HDAC3 Monoclonal Antibody
Catalog Number:	N0326
Reactivity:	Human
Applications:	WB,IF/ICC,ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human HDAC3 (aa224-428) expressed in E. Coli.
Concentration:	
Gene Name:	HDAC3
Protein Name:	Histone deacetylase 3
Other Name:	HDAC3; Histone deacetylase 3; HD3; RPD3-2; SMAP45
Human Swiss Prot No.:	O15379
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	HDAC3 Monoclonal Antibody detects endogenous levels of HDAC3 protein.
Storage and Stability:	-20°C/1 year