

Antibodies Product Datasheet

Product Name:	HFE Monoclonal Antibody
Catalog Number:	N0329
Reactivity:	Human
Applications:	WB,IF/ICC,ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human HFE expressed in E. Coli.
Concentration:	
Gene Name:	HFE
Protein Name:	Hereditary hemochromatosis protein
Other Name:	HFE; HLAH; Hereditary hemochromatosis protein; HLA-H
Human Swiss Prot No.:	Q30201
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	HFE Monoclonal Antibody detects endogenous levels of HFE protein.
Storage and Stability:	-20°C/1 year