

Antibodies Product Datasheet

Product Name:	HGK Monoclonal Antibody
Catalog Number:	N0330
Reactivity:	Human
Applications:	WB,ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of HGK (aa400-500) expressed in E. Coli.
Concentration:	
Gene Name:	MAP4K4
Protein Name:	Mitogen-activated protein kinase kinase kinase 4
Other Name:	MAP4K4; HGK; KIAA0687; NIK; Mitogen-activated protein kinase kinase kinase 4; HPK/GCK-like kinase HGK; MAPK/ERK kinase kinase 4; MEK kinase kinase 4; MEKKK 4; Nck-interacting kinase
Human Swiss Prot No.:	O95819
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	HGK Monoclonal Antibody detects endogenous levels of HGK protein.
Storage and Stability:	-20°C/1 year