

Antibodies Product Datasheet

Product Name:	HSP60 Monoclonal Antibody
Catalog Number:	N0339
Reactivity:	Human, Mouse, Rat, Monkey
Applications:	WB, IHC-p, IF/ICCFCM, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/200 - 1/1000. Immunofluorescence: 1/200 - 1/1000. Flow cytometry: 1/200 - 1/400. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human HSP60 expressed in E. Coli.
Concentration:	
Gene Name:	HSPD1
Protein Name:	60 kDa heat shock protein mitochondrial
Other Name:	HSPD1; HSP60; 60 kDa heat shock protein; mitochondrial; 60 kDa chaperonin; Chaperonin 60; CPN60; Heat shock protein 60; HSP-60; Hsp60; HuCHA60; Mitochondrial matrix protein P1; P60 lymphocyte protein
Human Swiss Prot No.:	P10809
Formulation:	Ascitic fluid containing 0.03% sodium azide, 0.5% BSA, 50% glycerol.
Specificity:	HSP60 Monoclonal Antibody detects endogenous levels of HSP60 protein.
Storage and Stability:	-20°C/1 year