

Antibodies Product Datasheet

Product Name:	Human IgG Monoclonal Antibody
Catalog Number:	N0343
Reactivity:	Human
Applications:	WB,ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Human IgG was isolated from human sera and purified by chromatography.
Concentration:	1 mg/ml
Gene Name:	
Protein Name:	
Other Name:	
Human Swiss Prot No.:	
Formulation:	Purified antibody in PBS containing 0.03% sodium azide.
Specificity:	Human IgG Monoclonal Antibody detects Human IgG protein.
Storage and Stability:	-20°C/1 year