

Antibodies Product Datasheet

Product Name:	Human IgG-Fc Monoclonal Antibody
Catalog Number:	N0344
Reactivity:	Human
Applications:	WB,IF/ICC,ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified fusion protein with human IgG (Fc Specific) tag.
Concentration:	
Gene Name:	
Protein Name:	
Other Name:	
Human Swiss Prot No.:	
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	Human IgG-Fc Monoclonal Antibody detects Human IgG-Fc protein.
Storage and Stability:	-20°C/1 year