

Antibodies Product Datasheet

Product Name:	IGFBP2 Monoclonal Antibody
Catalog Number:	N0354
Reactivity:	Human
Applications:	WB,IHC-p,IF(paraffin section),ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of protein IGFBP2 (aa180-328) expressed in E. Coli.
Concentration:	
Gene Name:	IGFBP2
Protein Name:	Insulin-like growth factor-binding protein 2
Other Name:	IGFBP2; BP2; IBP2; Insulin-like growth factor-binding protein 2; IBP-2; IGF-binding protein 2; IGFBP-2
Human Swiss Prot No.:	P18065
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	IGFBP2 Monoclonal Antibody detects endogenous levels of IGFBP2 protein.
Storage and Stability:	-20°C/1 year