

Antibodies Product Datasheet

Product Name:	IGF-IR Monoclonal Antibody
Catalog Number:	N0355
Reactivity:	Human
Applications:	WB,IHC-p,IF(paraffin section),ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of IGF-IR expressed in E. Coli.
Concentration:	
Gene Name:	IGF1R
Protein Name:	Insulin-like growth factor 1 receptor
Other Name:	IGF1R; Insulin-like growth factor 1 receptor; Insulin-like growth factor I receptor; IGF-I receptor; CD antigen CD221
Human Swiss Prot No.:	P08069
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	IGF-IR Monoclonal Antibody detects endogenous levels of IGF-IR protein.
Storage and Stability:	-20°C/1 year