

Antibodies Product Datasheet

Product Name:	IKK α Monoclonal Antibody
Catalog Number:	N0360
Reactivity:	Human
Applications:	WB,IF/ICC,FCM,ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. Flow cytometry: 1/200 - 1/400. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human IKK α expressed in E. Coli.
Concentration:	
Gene Name:	CHUK
Protein Name:	Inhibitor of nuclear factor kappa-B kinase subunit alpha
Other Name:	CHUK; IKKA; TCF16; Inhibitor of nuclear factor kappa-B kinase subunit alpha; I-kappa-B kinase alpha; IKK-A; IKK-alpha; I κ BKA; I κ B kinase; Conserved helix-loop-helix ubiquitous kinase; I-kappa-B kinase 1; IKK1; Nuclear factor NF-kappa-B
Human Swiss Prot No.:	O15111
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	IKK α Monoclonal Antibody detects endogenous levels of IKK α protein.
Storage and Stability:	-20°C/1 year