

Antibodies Product Datasheet

Product Name:	IKK β Monoclonal Antibody
Catalog Number:	N0362
Reactivity:	Human
Applications:	WB,IHC-p,IF(paraffin section),ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of IKK β expressed in E. Coli.
Concentration:	
Gene Name:	IKBKB
Protein Name:	Inhibitor of nuclear factor kappa-B kinase subunit beta
Other Name:	IKBKB; IKKB; Inhibitor of nuclear factor kappa-B kinase subunit beta; I-kappa-B-kinase beta; IKK-B; IKK-beta; IkbKB; I-kappa-B kinase 2; IKK2; Nuclear factor NF-kappa-B inhibitor kinase beta; NFKBIKB
Human Swiss Prot No.:	O14920
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	IKK β Monoclonal Antibody detects endogenous levels of IKK β protein.
Storage and Stability:	-20°C/1 year