

Antibodies Product Datasheet

Product Name:	IKK ϵ Monoclonal Antibody
Catalog Number:	N0363
Reactivity:	Human
Applications:	WB,ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of IKK ϵ (aa1-257) expressed in E. Coli.
Concentration:	
Gene Name:	IKBKE
Protein Name:	Inhibitor of nuclear factor kappa-B kinase subunit epsilon
Other Name:	IKBKE; IKKE; IKKI; KIAA0151; Inhibitor of nuclear factor kappa-B kinase subunit epsilon; I-kappa-B kinase epsilon; IKK-E; IKK-epsilon; IkbKE; Inducible I kappa-B kinase; IKK-i
Human Swiss Prot No.:	Q14164
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	IKK ϵ Monoclonal Antibody detects endogenous levels of IKK ϵ protein.
Storage and Stability:	-20°C/1 year