

Antibodies Product Datasheet

Product Name:	IL-11 Monoclonal Antibody
Catalog Number:	N0365
Reactivity:	Human
Applications:	WB,ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human IL-11 expressed in E. Coli.
Concentration:	1 mg/ml
Gene Name:	IL11
Protein Name:	Interleukin-11
Other Name:	IL11; Interleukin-11; IL-11; Adipogenesis inhibitory factor; AGIF; Oprelvekin
Human Swiss Prot No.:	P20809
Formulation:	Purified antibody in PBS containing 0.03% sodium azide.
Specificity:	IL-11 Monoclonal Antibody detects endogenous levels of IL-11 protein.
Storage and Stability:	-20°C/1 year