

Antibodies Product Datasheet

Product Name:	IL-34 Monoclonal Antibody
Catalog Number:	N0370
Reactivity:	Human
Applications:	WB,IF/ICC,FCM,ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. Flow cytometry: 1/200 - 1/400. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human IL-34 expressed in E. Coli.
Concentration:	
Gene Name:	IL34
Protein Name:	Interleukin-34
Other Name:	IL34; C16orf77; Interleukin-34; IL-34
Human Swiss Prot No.:	Q6ZMJ4
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	IL-34 Monoclonal Antibody detects endogenous levels of IL-34 protein.
Storage and Stability:	-20°C/1 year