

Antibodies Product Datasheet

Product Name:	IL-6 Monoclonal Antibody
Catalog Number:	N0371
Reactivity:	Human
Applications:	WB,ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human IL-6 expressed in E. Coli.
Concentration:	
Gene Name:	IL6
Protein Name:	Interleukin-6
Other Name:	IL6; IFNB2; Interleukin-6; IL-6; B-cell stimulatory factor 2; BSF-2; CTL differentiation factor; CDF; Hybridoma growth factor; Interferon beta-2; IFN-beta-2
Human Swiss Prot No.:	P05231
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	IL-6 Monoclonal Antibody detects endogenous levels of IL-6 protein.
Storage and Stability:	-20°C/1 year