

Antibodies Product Datasheet

Product Name:	IMP-3 Monoclonal Antibody
Catalog Number:	N0373
Reactivity:	Human
Applications:	WB,IHC-p,IF(paraffin section),ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human IMP-3 expressed in E. Coli.
Concentration:	
Gene Name:	IGF2BP3
Protein Name:	Insulin-like growth factor 2 mRNA-binding protein 3
Other Name:	IGF2BP3; IMP3; KOC1; VICKZ3; Insulin-like growth factor 2 mRNA-binding protein 3; IGF2 mRNA-binding protein 3; IMP-3; IGF-II mRNA-binding protein 3; KH domain-containing protein overexpressed in cancer; hKOC; VICKZ family member 3
Human Swiss Prot No.:	O00425
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	IMP-3 Monoclonal Antibody detects endogenous levels of IMP-3 protein.
Storage and Stability:	-20°C/1 year